



INTRAORAL CAMERA with autofocus and arm attachment

Ref. 100336



- ✓ 21.5-inch LCD touchscreen.
- ✓ High-definition (4K) camera.
- ✓ Multifunction buttons on the camera.
- ✓ Autofocus.
- ✓ Image and video recording to USB memory.
- ✓ QR code for downloading images via Wi-Fi.
- ✓ Bar mount.
- ✓ Android operating system.
- ✓ Display of 1 or 4 images on screen.
- ✓ Allows annotation on images in 2 colors.




TECHNICAL CHARACTERISTICS

- ✓ Camera: 4K
- ✓ Focus: 10–69 mm
- ✓ Monitor: 21.5 inches
- ✓ Height: 28 cm
- ✓ Width: 49 cm
- ✓ Depth: 4 cm
- ✓ Lens: 1/2.8" CMOS
- ✓ Viewing Angle: 105°
- ✓ Power Supply: 230 Vac, 12 Vdc
- ✓ Illumination: 8 LEDs
- ✓ Monitor Weight: 2525 g
- ✓ Camera Weight: 100 g
- ✓ Total Weight with Packaging: 3525 g



✓ Visualización en detalle.



✓ Visualización completa.



PACKAGE CONTENTS

- ✓ Bar mount (for the dental unit arm).
- ✓ Camera connection cable.
- ✓ Intraoral camera.
- ✓ Monitor.
- ✓ Camera mount.
- ✓ Camera covers.
- ✓ USB flash drive.
- ✓ 12V, 4A power adapter.



MESTRA

Talleres Mestraitua SL
 (+34) 944 530 388
 info@mestra.es - www.mestra.es
 48150 SONDIKA (Spain)



Management
 System
 ISO 9001:2015

www.tuv.com
 ID 1100018003