



MESTRA®

CHORRO DE ARENA AUTOMATICO
AUTOMATIC SANDBLASTER
 SABLEAUSE AUTOMATIQUE
 Ref. 080218

Nº	Ref.	Concepto	Description	Nº	Ref.	Concepto	Description
1	080240-14	Manorregulador	Pressure control	21	080270-31	Porta boquillas	Nozzle holder
2	080235-18	Temporizador	Timer	22	080235-08	Flexo pistola	Gun arm
3	080422-17	Interruptor luz	Switch	23	080010-05	Pasa cables	Cable bushing
4	080220-11	Manómetro	Pressure gauge	24	080235-27	Eje	Shaft
5	080218-07	Rejilla ventilación	Ventilation grill	25	030350-03	Espárrago	Stud
6	080218-05	Cristal	Glass	26	080550-15	Motor	Motor
7	080508-15	Barra de LED	LED bar	27	080218-10	Travesaño	Crossbar
8	080218-11	Sufrideras cesta	Basket rubbers	28	080218-09	Larguero izquierdo	Crossbar left
9	080218-04	Cesta	Basket	29	080218-08	Larguero derecho	Crossbar right
10	080253-10	Electroválvula 24V	Electrovalve 24V	30	080523-07	Tomillo punta eje	Shaft tip screw
11	080242-24	Filtro electroválvula	Electrovalve filter	31	080523-08	Tuerca eje	Shaft nut
12	080614-20	Fuente alimentación	Power source	32	080235-24	V-ring	V-ring
13	000101	Tubo Ø4x6	Tube Ø4x6	33	080515-08	Casquillo con fieltro	Sleeve with felt
14	000312	Racor codo plástico 1/8 Ø6	1/8 elbow Ø6 plastic connector	34	080218-15	Casquillo motor	Motor sleeve
15	080235-11	Racor prolongación inyector	Prolongation injector connector	35	080235-30	Junta tórica casquillo	Motor sleeve O-ring
16	080235-10	Inyector	Injector	36	080235-37	Bloqueo motor	Motor blockage
17	080235-09	Cuerpo pistola	Gun body	37	080218-06	Tapa cámara estancia	Sealed chamber lid
18	080270-48	Casquillo pistola	Guns sleeve	38	080228-03	Tapa inferior	Lower lid
19	080235-13	Manguera aspiración	Aspiration hose	38	000200	Cable	Cable
20	080270-26	Boquilla tungsteno	Tungsten nozzle	39	080218-23	Protector cristal (10 u)	Glass protector (10 u)