



HYGIENIC STORAGE CABINET

Ref. 100505

Continuous protection for instruments and work materials in clinics and laboratories.

In the day-to-day operations of a dental clinic, managing instruments after processing is a key link in maintaining good hygiene practices. Leaving tools exposed on work surfaces or in conventional drawers increases the risk of environmental contamination from aerosols and handling.

This storage chamber optimizes workflow in the dental surgery, offering a protected and isolated environment for instruments. As an added benefit, it builds maximum patient confidence, as they can visually perceive that the dentist implements additional barriers to maintain hygiene at all times.

Its compact and functional design -with dimensions similar to a household microwave- allows it to integrate perfectly into the surgery cabinetry or the clinic's organization area.



Key Features:

- ✓ **ULTRAVIOLET (UV) LIGHT TECHNOLOGY:** Acts continuously on surfaces, helping to passively keep the internal equipment free from environmental bacterial load.
- ✓ **DUAL TECHNOLOGY (UV + HEAT):** Combines the action of light with a constant thermal treatment (between 60 and 75°C) for an optimal preservation environment.
- ✓ **AUTOMATED CYCLE CONTROL:** Digital panel with a programmable timer and AUTO function (an optimized 60-minute cycle with an initial UV action followed by thermal drying). It allows you to monitor the process status at a glance.

- ✓ **INTEGRATED DRYING FUNCTION:** Efficient moisture management to prevent instruments from remaining in a humid environment.
- ✓ **PROFESSIONAL CONSTRUCTION:** High-resistance, easy-to-clean stainless steel interior and finishes, designed to maximize internal light reflection (mirror effect).

IMPORTANT NOTE: This equipment is intended exclusively for the organized storage and environmental protection of work instruments. It is not a medical device and under no circumstances replaces regulatory sterilization methods (autoclaving or other systems approved by current regulations).

Specifications:

Dimensions: 440 x 365 x 405 mm	Chamber material: Stainless steel
Weight: 6 kg	Electrical power: 460 W
Interior volume: 38 L	Power supply: 230 V, 50 Hz



MESTRA

Talleres Mestraitua SL
(+34) 944 530 388
info@mestra.es - www.mestra.es
48150 SONDIKA (Spain)



Sistema de
Gestión
ISO 9001:2015

www.tuv.com
ID 1100018003