



MESTRA

CLINIC Catalogue


**2
0
2
6**



www.mestra.es/en


Index

Instruments	03
Cleaning and disinfection	17
Machines	28
Furniture	33
Lightning	36
Alfabetical index	38



Spatulas

		Total length	Blade	Weight
1. Ref. 060001	Alginate spatula. Stainless steel (white handle)	205 mm	110 mm	0,061 kg
2. Ref. 060002	Plaste Ref. alginate spatula (Plastic) 12 u.	208 mm	90 mm	0,024 kg
3. Ref. 060003	Alginate spatula. Stainless steel (white handle)	190 mm	90 mm	0,037 kg
4. Ref. 060000	Alginate spatula	230 mm	110 mm	0,075 kg



Magnetic bar (instrument holder)

Ref. 050362

Easy to fit and may be secured by bolts to almost any surface.
(Spatulas not included).



Alginate mixing bowl

Ref. 060120

Diameter	110 mm	Capacity	400 cm ³
Height	78 mm	Weight	0,100 kg



Dappen glasses

Ref. 070161 Ref. 070162 Ref. 070163 Ref. 070164

Reference	Color
070161	Amber (12 u.)
070162	Blue (12 u.)
070163	Transparent (12 u.)
070164	Green (12 u.)

Ø 33 x (H) 35 mm | 0,042 kg (1 u.)
(L) 160 x (W) 110 x (H) 40 mm | 0,63 kg (12 u.)



Dappen glass support

Ref. 070165


Completely manufactured in polished stainless steel with transparent silicone legs to avoid sliding. It is very useful and convenient to organize the worktable. It has 6 different cavities for glasses. (Glasses not included).



Plastic articulator

(for dental correction models) 50 u.

Ref. 010110



Plastic articulator

for maxillary orthopedics)

Ref. 010130




Base former with magnet

Ref. 050110=SP Ref. 050120=IP
Ref. 050130=SG Ref. 050140=IG

INSTRUCTIONS FOR USE

1. Fit the base former with the metal plate.
2. Empty the impression as usual.
3. Insulate the model and fit the magnet on the metal plate.
4. Make the counterform that will remain fixed in the articulator pouring coating plaster on the model; let it polymerize. Then separate the model carefully from the counterform (which can be reused for another model).



Model formers (for orthodontics training)

S = Upper	I = Lower	S = Upper	I = Lower
Ref. 050060	Ref. 050062	Ref. 050064	Ref. 050066
0,081 kg	0,081 kg	0,073 kg	0,071 kg




Modelling carvers

Ref. 070320 (Set 6 units)


Ref. 070330 (Set 5 units)



Dam perforator




Clamp tongs



Clamp holder: sterilizable

Ref. 070110 Made of silicone. Sterilizable.





Plier holding support

Ref. 050364

A simple and practical tabletop component which helps to keep the work area tidy.

Cement spatulas


- ① **Ref. 060240** 178 mm 0,02 kg
- ① **Ref. 060242** 175 mm 0,02 kg
- ② **Ref. 060244** 175 mm 0,02 kg
- ③ **Ref. 060246** 178 mm 0,02 kg
- ④ **Ref. 060248** 175 mm 0,02 kg
- ⑤ **Ref. 060250** 195 mm 0,02 kg




Stainless steel holder. For instruments

Ref. 060235

Completely made of stainless steel.
Can store up to up to 7 working tools. (Spatulas not included)



Bending pliers




Plier Mathieu

Ref. 070356



Plier Weingart

Ref. 070354



Standard wire folding pliers (for wire cutting)
Ref. 070339





Dirty cotton container
Ref. 070155 (stainless steel)
Ø 60-51 x (H) 77 mm | 0,300 kg

Clean cotton container
Ref. 070160 (stainless steel)


Crown removing hammer

Ref. 070035-04
Long tip
Ref. 070035-05
Short tip



Amalgam applier

Ref. 070085	(L) 145 x (W) 26 x Ø 14 mm	0,031 kg
Ref. 070090	(L) 133 x (W) 30 x Ø 14 mm	0,027 kg




Tongs

Ref. 070380 = G	(L) 300 x (W) 75 x (H) 25 mm	0,125 kg
Ref. 070390 = M	(L) 230 x (W) 66 x (H) 17 mm	0,059 kg
Ref. 070400 = P	(L) 200 x (W) 72 x (H) 20 mm	0,059 kg



Mixing mortar

Ref. 060140
Made of stainless steel.



Matrix retainer (universal)

Ref. 070409



Cotton roll dispenser

Made of plastic sterilizable
at 135° C

Ref. 070061	(L) 120 x (W) 97 x (H) 55 mm 0,2 kg
(L) 120 x (W) 105 x (H) 60 mm 0,2 kg	
Ref. 070062	(L) 108 x (W) 100 x (H) 52 mm 0,3 kg
(L) 120 x (W) 105 x (H) 60 mm 0,3 kg	

Liquid doser

Ref. 070650

400 cc capacity bottle, ideal for dosing all types of liquids easily and accurately. It can be handled using just one hand.

Two scales: one in millimetres and the other in doses for preparing the alginate

Ø 77 x (H) 200 mm	0,061 kg
(L) 205 x (W) 80 x (H) 80 mm	0,102 kg



Glass dispenser

Ref. 070491	(L) 320 x (W) 100 x (H) 113 mm	0,354 kg
(L) 340 x (W) 135 x (H) 122 mm	0,490 kg	
Ref. 070490	(L) 358 x (W) 85 x (H) 85 mm	0,777 kg
(L) 373 x (W) 95 x (H) 95 mm	0,860 kg	



Bur holder (plastic)

Ref. 070142

Ø 90 x (H) 66 mm | 0,07 kg

(L) 91 x (W) 91 x (H) 71 mm | 0,09 kg



Bur holder (clinic)

Ref. 070150


Burs held by means of a magnet.
Made of (ferritic) stainless steel
with plastic casing.

(L) 115 x (W) 45 x (H) 45 mm | 0,18 kg



Ref. 070151 = FG

Ref. 070152 = CA



Bur holder (aluminium)

Ref. 070151 = FG | (L) 60 x (W) 25 x (H) 29 mm | 0,040 kg

Ref. 070152 = CA | (L) 40 x (W) 25 x (H) 27 mm | 0,046 kg

Magnetic bur holder

Ref. 070115

Very practical with an attractive design.
Its small size means it can be placed
anywhere.

Supplied in a batch of three products.
Colours: red, green
and blue. Ø 62 mm

Ø 62 x (H) 49 mm (1 u.) | 0,182 kg

(L) 200 x (W) 40 x (H) 60 mm - Paking (3 u.)
0,583 kg



Alcohol lamp

Ref. 070065

(stainless steel)

Ø 76 x (H) 54 mm | 0,100 kg

Lamp wicks:
Ref. 070065-03
Ref. 070065-10
Ref. 070065-11



Ref. 070070

(glass)

Ø 60 x (H) 83 mm | 0,192 kg

(L) 200 x (W) 40 x (H) 60 mm - Paking (3 u.)
0,583 kg

Descontaminadores de fresas

Ref. 070530

Ref. 070531



Descontaminadores de fresas

Ref. 070530

Ref. 070531



Only the base, made of temperature
resistant plastic. Capacity for 36 pieces.

Made for: FC, CA, FC+CA, laboratory.



Bur holder (plastic, sterilizable)

Ref. 070120 | (L) 75 x (W) 60 x (H) 35 mm | 0,060 kg

Ref. 070125 | (L) 75 x (W) 60 x (H) 35 mm | 0,060 kg

Ref. 070130 | (L) 75 x (W) 60 x (H) 35 mm | 0,060 kg


Ref. 070135 | (L) 75 x (W) 60 x (H) 55 mm | 0,072 kg

Ref. 070120 = FG

Ref. 070125 = CA

Ref. 070130 = FG+CA

Ref. 070135 = Laboratory




(L) 125 x (W) 5 x (H) 24 mm (1 u.) | 0,012 kg
(L) 130 x (W) 450 x (H) 27 mm Packing (4 u.) | 0,060 kg

Dental film clip (4 u.)

Made of plastic

Ref. 070194

Fully permeable to X-rays.
Enduring, reliable and easy to use.




(L) 38 x (W) 12 x (H) 8 mm (1 u.) | 0,002 kg
(L) 70 x (W) 50 x (H) 16 mm Packing (12 u.) | 0,017 kg

Dental film clip (12 u.)

Stainless steel

Ref. 070190



(L) 205 x (W) 65 x (H) 22 mm | 0,065 kg

Holder for 10 dental film clips

Ref. 070192



It simplifies the operation and handling of up to 10 dental films. Completely made in stainless steel.

X-ray viewers (flat)

X-ray viewers with a high brightness (3500 Lux) cold-light screen. This technology also features a low power consumption and a practically unlimited duration (10000 hours minimum).

Model Ref. 100132 has a large SRA3-size LED-backlit screen on which large-format films can be viewed. It also has light intensity regulation and two types of operation: automatic, when the film approaches the device turns on without pressing the button; switched on by pressing the button.


	Vision aerea	Outer area		
Ref. 100131 DIN A4	300 x 210 mm	360 x 240 x 20 mm	Wall and tabletop	1,700 kg
Ref. 100132 SRA3	430 x 355 mm	510 x 495 x 23 mm	Wall and tabletop	3,076 kg
Ref. 100134 DIN A2	580 x 400 mm	700 x 520 x 33 mm	Wall	4,000 kg



(L) 410 x (W) 200 x (H) 88 mm | 4,032 kg
(L) 515 x (W) 350 x (H) 290 mm | 4,500 kg

X-ray viewer **Ref. 070460**

Supplied with 30° inclined tripod. Can also be placed horizontally or screwed to the wall. Includes a system of springs by which the cover can hold the plates. Size 400 x 200 mm.



Dose measuring spoon **Ref. 100270**

PROVEN: Save 20-40% on plaster.

kg 0,042 kg



Ref. 100041

0.1 g resolution, maximum load of 3000 g, and "piece counting" function, very useful to count small elements such as pins, etc.

Ref. 100040

0.05 g resolution in the weight rank of 0-1000 g, and 0.1 g in the rank of 1000-2000 g.

Ref. 100042

Precise resolution of 0.1 g and high weighing capacity (1000 g).

Precision scale


Ref. 100040	(L) 150 x (W) 110 x (H) 36 mm	0,286 kg
Ref. 100041	(L) 248 x (W) 248 x (H) 100 mm	0,774 kg
Ref. 100042	(L) 300 x (W) 190 x (H) 55 mm	0,446 kg



Hand magnifying glass with light

Ref. 100092

Supplied with 4 set of lenses:
1.5x, 2.0x, 2.5x and 3.5x magnification



Magnifying glass with base


Ref. 100098

(L) 315 x (W) 154 x (H) 45 mm | 0,427 kg

Hand magnifying glass with light

Ref. 100099



(L) 70 x (W) 50 x (H) 30 mm | 0,06 kg




Holding clamps

Ref. 050554

(L) 120 x (W) 75 x (H) 68 mm
0,42 kg



Medidores de coronas y calibre



Microscope with base

Ref. 100097

5 - 15 magnification.
180 mm working distance. mm.

(L) 350 x (W) 275 x (H) 370 mm | 5 kg

Accessories (included):

- ✓ 2 slide holders (stainless steel)
- ✓ 3 glasses



Manual developer

Made of injected plastic with sun filter

Ref. 070170

(L) 520 x (W) 300 x (H) 360 mm | 2,8 kg



Manual developer (With negatoscope)

Ref. 070180

(L) 520 x (W) 300 x (H) 360 mm | 4,7 kg

Magnifying glasses

Adjustable system that adapts to the head and keeps the hands free to work. It has Ultra clear lens of minimal light absorption. It can be worn over glasses.

There are three different models available:

	Ref. 100079 (L) 285 x (W) 225 x (H) 89 mm	0,55 kg
	Ref. 100080 (L) 205 x (W) 185 x (H) 70 mm	0,25 kg
	Ref. 100090 (L) 175 x (W) 175 x (H) 70 mm	0,19 kg



Binocular magnifying glass

Ref. 100096 2,5 magnification.
420 mm working distance.




Galilean dental loupes **Ref. 100093**

Without flip-up mechanism, so they are lighter and the depth of field and field of vision is greater and wider than with the flip-up system. Four different types of loupes: SS, S, R, L.



LED light for galilean dental loupes **Ref. 100094**

- ✓ Snap on type headlight + yellow color filter.
- ✓ Switching adapter.
- ✓ Cordeless headlight controller (x2).
- ✓ Charging cable.



UV-ray safety screen

Ref. 070515



Safety glasses

Ref. 100011-12-13		(L) 160 x (W) 54 x (H) 56 mm	0,05 kg
Ref. 100014		(L) 180 x (W) 90 x (H) 65 mm	0,07 kg



Face protections

Ref. 100021 (10 u.)	Ref. 100022 (10 u.)
 (L) 285 x (W) 185 x (H) 35 mm 0,23 kg	 (L) 308 x (W) 215 x (H) 35 mm 0,28 kg



Magnifying glasses - Bioart

Magnification 2,5x Magnification 3,5x

Ref. 100076 **Ref. 100077**

Lightweight, handy and with excellent panoramic and clear vision. The magnifier is essential for clinical and laboratory procedures. Manufactured with CR-39 resin, higher resistance against abrasion and excellent and excellent transparency.

Made of plastic
sterilizable at 135° C

Ref. 070054



Ref. 070051

Ref. 070053

Cheek retractor

Ref. 070051		(L) 155 x (W) 10 x (H) 25 mm	Packing (2 u.)	0,324 kg
Ref. 070053		(L) 155 x (W) 10 x (H) 25 mm	Packing (2 u.)	0,355 kg
Ref. 070054		(L) 115 x (W) 70 x (H) 20 mm	Packing (2 u.)	0,234 kg

Lip retractor **Ref. 070055**

Made of plastic sterilizable at 135°

(L) 120 x (W) 53 x (H) 18 mm	(1 u.)	0,008 kg
(L) 145 x (W) 60 x (H) 25 mm - Paking	(2 u.)	0,028 kg



Lip retractor **Ref. 070056**

Made of plastic sterilizable at 135° C.

(L) 115 x (W) 42 x (H) 18 mm	(1 u.)	0,008 kg
(L) 140 x (W) 50 x (H) 25 mm - Paking	(2 u.)	0,026 kg



(L) 150 x (W) 45 x (H) 20 mm	(1 u.)	0,019 kg
(L) 170 x (W) 55 x (H) 30 mm - Paking	(2 u.)	0,038 kg




Ref. 070057 A	(L) 166 x (W) 55 x (H) 20 mm	Packing (1 u.)	0,019 kg
Ref. 070057 B	(L) 145 x (W) 45 x (H) 16 mm	Packing (1 u.)	0,012 kg
Ref. 070057 A	(L) 190 x (W) 60 x (H) 35 mm	Packing (2 u.)	0,055 kg
Ref. 070057 B	(L) 175 x (W) 50 x (H) 50 mm	Packing (2 u.)	0,044 kg

Lip retractor (2 units) double use

Ref. 070052

Made of plastic sterilizable at 135° C



Made of plastic sterilizable at 135° C

Lip retractor

Ref. 070058 A

Ref. 070058 B

Ref. 070058 C

Ref. 070058 A	(L) 126 x (W) 50 x (H) 53 mm	Packing (1 u.)	0,010 kg
Ref. 070058 B	(L) 105 x (W) 48 x (H) 48 mm	Packing (1 u.)	0,009 kg
Ref. 070058 C	(L) 98 x (W) 42 x (H) 40 mm	Packing (1 u.)	0,007 kg
Ref. 070058	(L) 145 x (W) 60 x (H) 55 mm	Packing (2 u.)	0,034 kg

Lip retractor

Ref. 070057A **Ref. 070057B**

Made of plastic sterilizable at 135° C



Napkin holder (10 u.)

Ref. 070060

(L) 520 mm (u.) | 0,014 kg (u.)

(L) 100 x (W) 120 x (H) 30 mm (aprox.)	0,161 kg (10 u.)
--	------------------



Drill cleaning tongs

Ref. 070520

Ideal for immersion of drills and other small items in disinfectant baths, ultrasonic equipment and autoclaves.

Made of stainless steel.

(L) 156 x (W) 45 x (H) 40 mm | 0,021 kg




Fox plane (plastic)

Ref. 070510

Fabricado en plástico
esterilizable a 135 C°.

(L) 180 x (W) 130 x (H) 6 mm | 0,038 kg




Cofferdam arc holder (young)

Ref. 070045

Of special heat-resistant plastics allowing it to be sterilized, unbreakable and radio transparent.

Ref. 070050 Metal.





Micro applicators (400 u. per box)

Ref. 070667 Ref. 070669 Ref. 070671

NEW



Disposable brushes for prophylaxis (100 u)


Ref. 070702 Hard Ref. 070704 Soft

NEW

Disposable cups for prophylaxis (100 u)

Ref. 070706 Hard Ref. 070708 Soft

NEW



Mirror handle

Ref. 070010
Stainless steel.

Ref. 070020
Anodized aluminium.
Colour: Silver.

Ref. 070030
Stainless steel, securely clamping
the points as the clamp is divided
into three sections.



Assorted toy bag

Ref. 160000

Assortment of 25 units.
Ball size: Height 550 mm x Width 500 mm

NEW

Mouth mirrors (12 units)

Top quality mirrors meeting all the required standards.
They all are plane with stainless steel chassis. They are
available in silver (normal) and rhodium.

The image is more clear and with no distortions.
There are four different sizes: Plane N°3, N°4, N°5 and N°6.
These also have two kinds of thread: cone socket
and simple stem.

Cone socket

Plane	Normal	Rhodium
N° 4	070551	070560
N° 5	070552	070561
N° 6	070553	

Simple stem

Plane	Normal	Rhodium
N° 4	070556	070565
N° 5	070557	070566



Packing (12 u.) (L) 84 x (W) 52 x (H) 28 mm | 0,080 kg


Mixing tips

Supplied in bags of 50 pcs.

These mixing tips are compatible with the main injection systems in the market.

REFERENCE	COLOR
Ref. 100172	Yellow (50 uds.)
Ref. 100176	Yellow Points (50 uds.)
Ref. 100170	Green (50 uds.)
Ref. 100178	Blue claro larga (50 uds.)
Ref. 100180	Blue claro corta (50 uds.)
Ref. 100184	Pink (50 uds.)
Ref. 100187	Brown (50 uds.)
Ref. 100187-01	Points 1,06 mm (50 uds.)
Ref. 100187-02	Points 2,00 mm (50 uds.)
Ref. 100185	Brown short (50 uds.)
Ref. 100186	Brown long (50 uds.)

(L) 160 x (W) 110 x (H) 45 mm | 0,12 kg



Silicone pistol
(Reinforced)

Ref. 100166

(L) 190 x (W) 120 x (H) 70 mm | 0,22 kg



Disposable trays (100 u.)

Ref. 100219 Complet kit

Hygienic and disposable trays. They save time and detergent because they cut out cleaning operations. They are adjusted to the specific dimensions of each patient by applying heat.

Transparent, so they can be used in implantology.

Available in 10 different shapes and sizes.

Content of the box:

10 units - UL - UM - US - P1
10 units - LL - LM - LS - P2
10 units - PB - PA

Replacement: 25 units of the ordered model

Ref. 100220-UL - Upper large
Ref. 100220-UM - Upper medium
Ref. 100220-US - Upper small
Ref. 100220-LL - Lower large
Ref. 100220-LM - Lower medium
Ref. 100220-LS - Lower small
Ref. 100220-PA - Partial A
Ref. 100220-PB - Partial B
Ref. 100220-P1 - Lower (kid)
Ref. 100220-P2 - Upper (kid)

(L) 380 x (W) 290 x (H) 130 mm Paking (100 u.) | 1,526 kg

(L) 200 x (W) 200 x (H) 80 mm Paking (25 u.) | 0,246 kg



Disposable trays (100 u.) with retention

Ref. 100222 Complet kit


Made of a very resistant plastic, and therefore rigid and non-deformable. Thanks to their anatomic design and to the different sizes and shapes available (10 options) they can be easily adaptable to any patient.

They feature inner retention to prevent unwanted movements of the alginate. Besides, their base has several cuts so that -if necessary- windows can be easily opened with a drill, thus avoiding the implant posts.

Content of the box:

10 units UL1 - UM3 - US5
10 units LL2 - LM4 - LS6

10 units Q7 - Q8
10 units UA9 - LA10



Replacement: 25 units of the ordered model

Ref. 100222-UL1 Upper large
Ref. 100222-UM3 Upper medium
Ref. 100222-US5 Upper small
Ref. 100222-LL2 Lower large
Ref. 100222-LM4 Lower medium
Ref. 100222-LS6 Lower small
Ref. 100222-Q7 Partial 7
Ref. 100222-Q8 Partial 8
Ref. 100222-UA9 Lower front
Ref. 100222-LA10 Upper front

(L) 380 x (W) 290 x (H) 130 mm | 1,526 kg

Impression trays (with retention) Stainless steel

Ref. 070450 Set of 6 trays


Ref. 070452 = L2 0,047 kg Ref. 070457 = U2 0,055 kg
 Ref. 070453 = L3 0,048 kg Ref. 070458 = U3 0,052 kg
 Ref. 070454 = L4 0,043 kg Ref. 070459 = U4 0,049 kg



Impression trays - Implants (10 u. of each reference)

Ref. 100222


BLUE	100222-CA	Partial	A
	100222-CB	Partial	B
	100222-SG	Upper	Large
	100222-IG	Lower	Large
	100222-SM	Upper	Medium
	100222-IM	Lower	Medium
	100222-SP	Upper	Small
	100222-IP	Lower	Small
	100222-SA	Upper	Anterior
	100222-IA	Lower	Anterior



Impression trays - Autoclavable (12 u. of each reference)

Ref. 100223

GREEN	100223-CA	Partial	A
	100223-CB	Partial	B
	100223-SG	Upper	Large
	100223-IG	Lower	Large
	100223-SM	Upper	Medium
	100223-IM	Lower	Medium
	100223-SP	Upper	Small
	100223-IP	Lower	Small
	100223-SA	Upper	Anterior
	100223-IA	Lower	Anterior



Fluoride foam trays (50 u. of each reference)

Ref. 100224

BLUE	100224-G	Double	Large
YELLOW	100224-M	Double	Medium
GREEN	100224-P	Double	Small



Transport boxes

Ref. 070605


Cord for boxes.
Supplied in 5 assorted colours
(25 units). Length: 1 meter

Ref. 070600

(25 u.) (L) 300 x (W) 180 x (H) 160 mm
0,74 kg (25 u.)



Inner dimensions



(25 u.) (L) 240 x (W) 160 x (H) 140 mm
0,54 kg (25 u.)




Orthodontic and denture boxes

Made of plastic. 25 units. Assorted colours

Ref. 070600

Inner dimensions




Boxes for orthodontic plates

Made of plastic. 25 units. Assorted colours

Ref. 070590

Inner dimensions




(100 u.) (L) 420 x (W) 380 x (H) 240 mm | 3,420 kg (100 u.)
(25 u.) (L) 300 x (W) 180 x (H) 160 mm | 0,742 kg (25 u.)



Ref. 070000-25 (25 u.)
Ref. 070001-25 (25 u.)

Inner dimensions



Ref. 070002-25 (25 u.) | **Ref. 070006-25** (25 u.)
Ref. 070007-25 (25 u.) | **Ref. 070008-25** (25 u.)
(100 u.) (L) 345 x (W) 305 x (H) 195 mm | 0,23 kg (100 u.)
(25 u.) (L) 240 x (W) 160 x (H) 140 mm | 0,23 kg (25 u.)



Denture boxes

Made of plastic


Ref. 070000 Blue (100 u.) **Ref. 070001** Pink (100 u.)

Boxes for orthodontic plates


Made of plastic in different colours (100 u.)

Ref. 070002 Blue **Ref. 070007** Light green **Ref. 070006** Green **Ref. 070008** Orange

Inner dimensions



Inner dimensions



Boxes for orthodontic plates

Made of plastic. 5 u. assorted

Ref. 070622

NEW


Boxes for orthodontic plates

Made of plastic. 5 u. assorted

Ref. 070616

NEW

Inner dimensions



Square boxes for orthodontic plates

Made of plastic. 5 u. assorted

Ref. 070618 White


Ref. 070619 Blue

Ref. 070620 Black

NEW

Boxes for orthodontic models

Ref. 070570 (100 u.)



(100 u.) (L) 555 x (W) 390 x (H) 335 mm | 12,330 kg (100 u.)



Container for models

Ref. 070003 horizontal 6 models

Ref. 070004 vertical 12 models

Ref. 070003 (L) 230 x (W) 240 x (H) 58 mm | 0,345 kg

Ref. 070004 (L) 290 x (W) 220 x (H) 75 mm | 0,425 kg

Transport boxes

Boxes with membrane

Made of plastic. Safety lock

Specially designed for the transport and delivery of high quality elements.

Thanks to its transparency, the content can be seen from every angle and perspective. Supplied in bags of 10 units.

The inside is cushioned by a high flexibility membrane, which gives you the impression that the content inside is hanging in the air.

Ref. 100262	(L) 55 x (W) 49 x (H) 25 mm	0,023 kg u.
10 u.	(L) 147 x (W) 56 x (H) 56 mm	0,261 kg
Ref. 100263	(L) 60 x (W) 100 x (H) 22 mm	0,045 kg u.
10 u.	(L) 125 x (W) 105 x (H) 135 mm	0,475 kg
Ref. 100264	(L) 100 x (W) 100 x (H) 45 mm	0,081 kg u.
10 u.	(L) 430 x (W) 110 x (H) 100 mm	0,911 kg



Boxes with membrane (10 unit bags) without safety lock

The interior is cushioned by a highly flexible membrane. Being transparent, the contents can be viewed from all angles and perspectives. Supplied in bags of 10 units.


Ref. 100259	(L) 23 x (W) 23 x (H) 25 mm	0,026 kg u.
Ref. 100260	(L) 50 x (W) 50 x (H) 25 mm	0,023 kg u.
Ref. 100261	(L) 75 x (W) 75 x (H) 40 mm	0,045 kg u.



Brackets boxes

Ref. 100263B (10 unit bags) Made of plastic

Ref. 100263 (10 u.) | (L) 60 x (W) 100 x (H) 22 mm | 0,045 kg u.



Boxes for models


Made of plastic

Small box

Ref. 070033 (50 u.)

Large box

Ref. 070034 (50 u.)



Ref. 070033	Packing (50 u.)	(L) 290 x (W) 225 x (H) 255 mm	1,121 kg
Ref. 070034	Packing (50 u.)	(L) 290 x (W) 225 x (H) 255 mm	0,930 kg



Transport and cleaning boxes

Ref. 100265 Made of plastic


(L) 115 x (W) 95 x (H) 46 mm	0,118 kg
(L) 120 x (W) 100 x (H) 50 mm	0,143 kg



Boxes for prosthesis (with drain) (4 units)

Ref. 100267 Made of plastic

Ref. 100267	(L) 85 x (W) 76 x (H) 80 mm	0,082 kg
4 units	(L) 165 x (W) 165 x (H) 80 mm	0,389 kg




Transport boxes

Made of plastic

Ref. 100255	(L) 29 x (W) 29 x (H) 23 mm	0,005 kg u.
50 u.	(L) 178 x (W) 150 x (H) 44 mm	0,289 kg- 50 u.
Ref. 100256	(L) 56 x (W) 44 x (H) 23 mm	0,011 kg u.
24 u.	(L) 280 x (W) 235 x (H) 23 mm	0,318 kg- 25 u.





Ref. 070444 Small	(L) 185 x (W) 145 x (H) 40 mm	0,241 kg
Ref. 070445 Large	(L) 290 x (W) 190 x (H) 40 mm	0,400 kg

Instrument tray


Made of anodized aluminium

Ref. 070444 Small

Ref. 070445 Large

SMALL


Ref. 070430	Silver
Ref. 070431	Gold
Ref. 070432	Red
Ref. 070433	Blue
Ref. 070434	Green



Ref. 070430-1-2-3-4 | (L) 185 x (W) 145 x (H) 40 mm | 0,250 kg

LARGE

Ref. 070439	Green
Ref. 070440	Silver
Ref. 070441	Gold
Ref. 070442	Red
Ref. 070443	Blue



Ref. 070439-0-1-2-3 | (L) 290 x (W) 185 x (H) 40 mm | 0,414 kg

Instrument tray

Made of anodized aluminium



Stainless steel box 18/8.
Anodized aluminium plate.

Ref. 070410 | (L) 210 x (W) 115 x (H) 45 mm | 1,127 kg



Ref. 070420 | (L) 85 mm | 0,017 kg

Box for endodontics

Ref. 070410

Ruler (for calibrating files)

Ref. 070420

Glazing furnace

Glazing furnace

Ref. 100252

What you've been waiting for

WHAT IS GLAZING AIR?

Glazing is the process of post-processing the sintered zirconia. Not only does it give a high aesthetic by forming a gloss layer on the prosthesis, but it also stabilizes the color of the prosthesis and increases the surface hardness to increase the resistance to tooth abrasion and fracture.

Zirconia prosthetics have become the trend, and more customers have to buy a porcelain furnace just for glazing.

The **GLAZING AIR** is furnace only for zirconia glazing, which boldly excludes vacuum functions not used for glazing, and is a affordable product that leaves only fast and precise temperature control technology, compact size and economic cost.

ADVANTAGES

AFFORDABLE FURNACE ONLY FOR GLAZING



Furnace only for zirconia glazing that excludes functions such as vacuum, and leaved absolutely necessary functions.

SPEEDY GLAZING



It reaches 1000 °C in less than 10 (average 5 minutes) minutes using a Kanthal® heater with good thermal conductivity.

PRECISE TEMPERATURE CONTROL



Precise temperature control provides the perfect outcome for any brand.

HIGH USER CONVENIENCE



- Both the membrane keyboard and the touch screen can be used simultaneously.
- Improves user convenience by including an auxiliary pedestal in the design (to place the result).

Rating power	Máx. 1,3 kW
Input power source	AC 220~240V 50/60Hz
Fuse	10 A
Max. temperature	1000 °C



Cleaning and disinfection

Autoclave 18-D, class B 18 L

Ref. 100365




- ✓ Class B 18-liter desktop autoclave.
- ✓ High quality internal components of Italian, Japanese and American brands.
 - ODE valves.
 - ULUAK vacuum pump.
 - ULKA water pump.
- ✓ Ultra quiet (noise level below 70 dB).
- ✓ Triple fractionated vacuum system with an independent steam generator, which guarantees fast and safe sterilization cycles.
- ✓ High precision control system for perfect sterilization results (EN13060).
- ✓ Self-diagnosis with error codes.
- ✓ Fully automatic with pre-programmed sterilization cycles.
- ✓ Stainless steel chamber and door, with electro-polished finish.
- ✓ Open type water tank, which allows easy filling and cleaning.
- ✓ LCD screen, which shows information about the cycle process: time, temperature, pressure, process alert.
- ✓ Different languages to choose.
- ✓ Print the sterilization report and display available in several languages.
- ✓ Includes USB and printer.
- ✓ Test: Vacuum test / Bowie & Dick test / Helix test.
- ✓ Pressure safety valve prevents overpressure in the chamber and steam generator.

Voltage	220-240 V 50/60 Hz.
Circuit breaker	F16 A / 400 V
Nominal power	1750 V
Sterilization temperatures	121 / 134 °C
Distilled water tank capacity	2.5 litres
Operating temperature (°C)	5-40
Relative operating humidity	Max. 80 %, not condensed
Max. level of noise	< 70 dB
Atm. operating pressure (kPa)	76-106

(L) 690 x (W) 575 x (H) 525 mm | 51 kg

icanCLAVE



Camera dimensions 247 / 350 mm



Filling tank.



Safety valve



Safety lock



Bacteriological filter

REGULATIONS

- ✓ 97/23/EEC Pressure Equipment Directive.
- ✓ Regulation (EU) 2017/745 on medical devices (MDR)
- ✓ EN13060 Small steam sterilizers.

- ✓ EN61326, EN 61000 Electromagnetic compatibility.
- ✓ EN61010-1, UL 61010-1 Safety requirements.
- ✓ EN61010-2-040, IEC 61010-2-040 Specific requirements for steam sterilizers used for treatment of medical materials.

- ✓ EN 13485:2012/AC:2012
- ✓ ISO 13485:2003

Official technical service ICANCLAVE in Spain

Cleaning and disinfection

Autoclave icanCLAVE 3F

Ref. 100367



icanCLAVE 3F is super-fast autoclave, especially for the sterilization of unwrapped instruments, which greatly improves the turnover rate of surgical instruments and the convenience of small operations as dentistry.

- ✓ **Printer port:** It is used for connecting an external printer print the cycle report automatically.
- ✓ **USB port:** it is used for connecting a USB disc, all sterilization data are automatically output in the USB disc and can be ready directly in any PC and stored electronically.
- ✓ **Internal Memory:** The last 9999 cycles can be automatically stored in the autoclave system, which can be printed out at any time.
- ✓ **Barcode label printer port:** It is used for connecting an optional barcode printer prints barcode labels that contain the information related to the sterilization cycle.
- ✓ **Suitcase for transport.**

CHARACTERISTICS:

- ✓ **Standard pre-programmed cycles:** Quick S, Standard B, Soft B, Prion 134°C, B&D Test, Helix Test, Vacuum Test.
- ✓ **Other cycles:** Descalc (cleaning steam generator and tube).
- ✓ Steam generator
- ✓ External water reservoir and waste water reservoir.
- ✓ Door system.
- ✓ Color touch screen.
- ✓ User friendly interface.
- ✓ Warning system.
- ✓ Atmospheric pressure identification system.
- ✓ **Safety system:**
 - Over pressure safety valve.
 - Over temperature switch.
 - Safety lock system preventing improper door lock.
 - Circuit breakers protecting the machine when it has short-circuited or over current.
 - Automatic warning for maintenance.




Touch screen.

Model Number	STE-3-F
Chamber dimension (mm)	192(W) x 45(H) x 285(D)
Chamber Volume	Aprox. 3L
Number of Trays	1
Material of Trays	Aluminum
Size of Trays (WxH)	178 x 285mm
Voltage. Freq.	220/110V, 50/60Hz
Power	2300 W

REGULATIONS:

- 93/42/EEC Directive on Medical Devices.
- Pressure Equipment Directive 2014/68/EU.
- EN 13060 Small Steam Sterilizers.
- EN 61326, EN 61000 Electromagnetic Compatibility.
- EN 61010-1, UL 61010-1 Safety requirements.
- EN 61010-2-040, IEC 61010-2-040.
- RoHS Directive 2011/65/EU.

icanCLAVE



Suitcase for transport.



Cleaning and disinfection



Microblasting unit (Laboratory)

Ref. 100122

They are essential when cleaning and degreasing microretentions, crowns, occlusal sandblasting, adjustment of interproximal contacts and a great many other precise sandblasting operations. **This model is for dental laboratory use; it cannot be used for intraoral applications.**

Tank for distilled water

Ref. 080620

A tank specially designed for the storage and supply of distilled water under pressure in the dental sector. Made of stainless steel, it features a pressure regulator and a double safety valve.

Before filling the tank, it must be properly cleaned and sterilized. Do not use detergents containing hydrochloric acid, hypochlorous acid, etc.

- ✓ Pressure: 5 bar
- ✓ Connection: Ø 6 mm.
- ✓ Made of stainless steel
- ✓ Pressure regulator and a double safety valve.

Capacity	18 liters
Weight	4,7 kg



Biological spores incubator

Simple and easy to use spore incubators. They are essential to carry out the necessary tests that confirm the proper performance of the sterilization devices. Two models are available: **Elite**, with digital display and adjustable temperature; and **Practic**, very functional and economical. Both have four cameras to be able to perform as many verifications at the same time.

The **Elite** model also has a screen that indicates the temperature at which the test is performed and the time remaining until it ends. Its electronic control also alerts of anomalies such as power outages to avoid incorrect results.

Both devices perform the test at 57,5°C, although the Elite model also allows operation at 35°C depending on the type of spores used.

35°C +- 2,5°C for Bacillus Subtilis
57,5°C +- 2,5°C for Geobacillus Stearothermophilus

Reference	Élite Ref. 100370	Practic Ref. 100375
Voltage	230 V 50/60 Hz	90-250 V 50/60 Hz
Consumption	25 W	10 W
Operating temperature	57,5 +- 2,5°C e 35 +- 2,5°C	57,5°C
Height	80 mm	90 mm
Width	105 mm	75 mm
Depth	165 mm	75 mm
Weight	0,750 kg	0,239 kg

Ref. 100370	(L) 210 x (W) 260 x (H) 175 mm	1,1 kg
Ref. 100375	(L) 103 x (W) 110 x (H) 95 mm	0,3 kg



Ref. 100370

Ref. 100375

Distiller/decalcifier

Ref. 100210

Simply fill the tank with tap water to get 4 litres of distilled water in around 5 hours. The cleaning of the tank is easy and quick. The water produced is very pure (between 2 and 10 mg/ litre of solid residue) and can be used in a large number of applications in the dental clinic and laboratory (autoclaves, creation of compounds, steam generating machines, etc.).

The appliance is supplied with a 4-litre PVC tank, a set of activated carbon filters and a can of citric acid to clean the tank.

IMPORTANT: The water supplied for the appliance is "DISTILLED". However, it cannot be used for syringes, vials, cleaning wounds and generally any surgical application that requires TRI-DISTILLED purity.

Ref. 100210-01

Distiller cleaner refill (500 g.)

Ref. 100210-02

Replacement filters for distiller (6 units)



Dimensions	280 x 250 x 365 mm
Weight	3,5 kg
Voltage	230 V 50/60 Hz
Power	580 W
Capacity	4 l.

(L) 285 x (W) 285 x (H) 400 mm | 5 kg



Cleaning and disinfection

Ultrasonic cleaner DA10

Ref. 100750

HIGH-PERFORMANCE ELECTRONICS

- ✓ Equipped with a state-of-the-art chipset and advanced algorithms.
- ✓ Automatically distinguishes between tooth and gum.
- ✓ Ensures a comfortable and safe cleaning experience.

HIGH QUALITY HANDPIECE WITH ILLUMINATION

- ✓ Removable design allows for high temperature and high pressure sterilisation.
- ✓ Tinted light for improved visibility during treatment.

INTELLIGENT POWER LEVEL FEEDBACK

- ✓ Automatically adjusts the power according to the hardness of the stone.
- ✓ Increased cleaning efficiency.



Weight	1,3 kg
Input	35VA
Output	3W-29W
Frequency	18kHz-45kHz
Voltage	230 V, 50 Hz

Ultrasonic cleaner Q5


Ref. 100755

- ✓ Power setting (100-800 Mn).
- ✓ **3 operating modes:** General cleaning, Periodontal treatment, Endodontic treatment.
- ✓ Automatic water supply from a bottle integrated in the unit.
- ✓ **Waterproof control panel.**
- ✓ **Large capacity water tank (800 ml)**
- ✓ **Includes IPX4 Wireless footswitch:** Eliminates tangled wires.

SINUSOIDAL WAVE TECHNOLOGY

Avoids the instantaneous impact of the traditional square wave, providing stable amplitude and regular, powerful vibrations. Efficiency and comfort for painless tooth cleaning.

PAINLESS TECHNOLOGY
Minimize patient sensitivity.



INCLUDED ACCESSORIES

- ✓ LED PH-2 handpieces x4
- ✓ Cordless footswitch x1
- ✓ Corded footswitch x1
- ✓ Power cord x1
- ✓ Water bottle x1
- ✓ Endodontic Wrench x1
- ✓ Cleaning Tips x16
- ✓ Waterproof Cover Pad x1
- ✓ Endodontic Files #25 x12

Weight	1,8 kg
Input	35VA
Output	3W-20W
Frequency	25kHz-35kHz
Voltage	230 V, 50 Hz

Ultrasonic cleaner Q6

Ref. 100760

6 functions in 1 device:

- ✓ Ultrasonic cleaning.
- ✓ Periodontal therapy.
- ✓ Root canal irrigation.
- ✓ Supragingival air polishing.
- ✓ Subgingival air polishing.
- ✓ Implant maintenance.

TOUCH PANEL:


Clear interface with ultra-sensitive sensors, dust and water resistant design. Works even with gloves, water drops or protective covers.

DUAL WATER TANK DESIGN:

Quickly switch between tanks to use pure water or medicated solutions with a single click.

TITANIUM ALLOY TIPS

Low hardness to protect teeth and gums.



Weight	1,5 kg
Input	30 VA ~ 48 VA
Output	3 W ~ 20 W
Frequency	25kHz-35kHz
Voltage	100 V-240 V~, 50/60 Hz

PAINLESS TECHNOLOGY
Minimize patient sensitivity.

Cleaning and disinfection

Basic curing lamp V200

Ref. 100700


The robust fibre optic cannon effectively directs the light to achieve homogeneous **polymerisation of the composites.**

Its outstanding value for money, making it an ideal choice for new clinics, or as an indispensable second unit, offering peace of mind in the event of unforeseen events.

FEATURES

- ✓ Constant light intensity.
- ✓ Up to 600 polymerisations with a single recharge.
- ✓ Comfortable and ergonomic design. Easy to use.
- ✓ Automatic switch-off in case of non-use.
- ✓ Radiation: between 1600 ~ 2000 W/cm²
- ✓ Emission limit 3.9 X 10J
- ✓ Polymerising wavelength: 420~480nm.

Dimensions	243 mm x Ø25 mm
Weight	390 g
Capacity	2200mAh
Voltage	3.7V
Battery	Rechargeable



Mid-range curing lamp DV 50


Ref. 100705

Designed to optimise every step of composite polymerisation.

The main innovation of the DV50 lies in the direct placement of the LED emitter. This results in a significantly higher power density, reaching up to 2600 mW/cm². As a result, curing takes place very quickly, minimising exposure time.

Ergonomic and practical. Its straight, streamlined shape allows effortless access to posterior areas of the arch, such as molars and third molars, eliminating the discomfort associated with other lamps.

Dimensions	217 mm x Ø 22,5 mm
Weight	117 g
Maximum intensity	2600mW/cm ²
Wavelength	385 nm-515 nm
Battery	Rechargeable 2200mAh
LED	Emitter of 10 W
Power adapter	≈ 100-240 V 50/60 Hz.
Luminous intensity measurement function	



High-end curing lamp VAFU

Ref. 100710

A device in the premium range of dental polymerisation. Its light power reaches 3200 mW/cm², which allows ultra-fast curing in just 1 or 2 seconds, drastically minimising the time the lamp remains in the mouth and optimising patient comfort.

exceptionally broad spectrum of light (385 nm to 515 nm). This ensures efficient polymerisation in record time of virtually all dental materials available on the market.

Advanced multi-LED head, consisting of four independent LEDs. Each emits at a different frequency, achieving an

Dimensions	33 x 27 x 206 mm
Power	3200 mW/cm ²
Wavelength	385nm-515nm
Battery	Rechargeable 2200 mAh.
Power adapter	≈ 100-240 V 50/60 Hz.
Luminous intensity measurement function	



INCLUDES:

- ✓ 1 main unit.
- ✓ 1 charger.
- ✓ 1 Batterie.
- ✓ 1 set of lenses (4 pcs.)
- ✓ 1 light protection shield.
- ✓ 1 cavity detector head.
- ✓ Goggles.


Cleaning and disinfection

Automatic bag sealer

Ref. 080435

- ✓ 7" touch screen with indications.
- ✓ Quick heating system for a perfect seal.
- ✓ The sealing temperature can be regulated by the user.
- ✓ Its integrated printer allows marking the obtained records, for a correct traceability, both in the clinic and in the laboratory.
- ✓ In accordance with standards ISO11607-2: 2006 and EN868-5-2009

Start Date: 20200512 Exp. Date: 20200611 Batch No.: 00000000 Operator: Higienista Department: Sección esteril Nota: Prog. envuelto ID: 5119C032



Sealing speed	10 m/minute
Temperature	60-220°C
Sealing pressure	86 N
Power	500 W
Dimensions	635 x 240 x 200 mm
Weight	17 kg

Ref. 080436
Optional roller



NEW

Bag sealer


Ref. 080433

Simple and nice looking bag sealer.

Despite its small size, it has a sealing length of 305 mm and a sealing width of 12 mm. These features are similar to those of sealers of much higher cost. It is ideal for the dental clinic.

Height	210 mm
Width	430 mm
Depth	335 mm
Weight	6.5 kg
Voltage	230 V 50/60 Hz
Power	100 W

(L) 500 x (W) 415 x (H) 245 mm | 8,7 kg



Barrier film

Ref. 080439

This adhesive film offers a practical and effective solution for dental environments, providing a layer of protection for instruments and materials.

Each roll contains 1200 pre-cut 10 x 15 cm adhesives. Blue colour.



Ref. 080447

Ref. 080446

Ref. 080445

Ref. 080444

Ref. 080443

Ref. 080441

Ref. 080442

Sterilization rolls

Ref. 080441	5 cm x 200 m
Ref. 080442	7.5 cm x 200 m
Ref. 080443	10 cm x 200 m
Ref. 080444	15 cm x 200 m

Ref. 080445	20 cm x 200 m
Ref. 080446	25 cm x 200 m
Ref. 080447	30 cm x 200 m

Class I medical device (MDR 2017/745)
Conforming to standards ISO 11607-1:201 and EN 868-5:2009

New bags 5 x 9 cm, 200 units per package
THE SMALLEST ON THE MARKET!

NEW



Self-sealing sterilisation pouch

Ref. 080448	5 x 24 cm
Ref. 080449	7 x 26 cm
Ref. 080450	9 x 26 cm
Ref. 080451	5 x 9 cm
Ref. 080452	14 x 28 cm

Ref. 080454	19 x 36 cm
Ref. 080456	25 x 40 cm
Ref. 080457	30 x 40 cm
Ref. 080458	30 x 45 cm
Ref. 080459	5.7 x 13 cm

Class I medical device (MDR 2017/745)
Conforming to standards ISO 11607-1:201 and EN 868-5:2009

Cleaning and disinfection



Ultrasonic cleaning device for splints

Ref. 100162

Ultrasonic cleaning device for dental elements, with a UV lamp that is used to disinfect during cleaning. Its 28 W of power ensures excellent results.

Its small dimensions make it suitable for home environments and travel. It also has a digital timer.




Ultrasonic bath 2,5 l.

Ref. 100140

Very comfortable and practical for cleaning small objects in the dental laboratory and Clinic. Its 2.5-litre capacity tank is the right size for holding the most commonly used dental laboratory instruments. It has a heating system and an electronic timer.

(L) 360 x (W) 275 x (H) 240 mm | 3,5 kg



Ultrasonic bath 1,4 l.

Ref. 100141

Essential machine in both dental lab and clinic. Real power of 60 W. that guarantees an excellent cleaning power. It has an electronic timer with 7 segments display.

(L) 430 x (W) 270 x (H) 270 mm | 2,5 kg



Ultrasonic bath (mini)

Ref. 100144



Ultrasonic bath

Ref. 100160 **Ref. 100160-10** (10 u.)
For denture cleaning.



Ultrasonic bath, 3 l.

Ref. 100142

For cleaning small objects in the laboratory and dental practice. With detachable cable, drain valve and hose. Includes plastic protective basket and drain hose.



Cleaning and disinfection



Calypso 4l. ultrasonic cleaner

Ref. 080422

The "Calypso" model is a last generation Ultrasonic Cleaner that includes the latest innovations: Sandwich type transducers, with 100 W real power transmitted to the tank; working frequency of 40 KHz., representing an optimum compromise between sonority and cleaning power; electronic time control run by a microprocessor, with digital display and warning acoustic signals; automatic heating system of the cleaning liquid; high density polyurethane body, a material resistant to chemical agents; tank of 4 litre capacity and waste pipe; etc. All these characteristics are complemented by an attractive line of careful design.

The machine is supplied with the following accessories:

- ✓ A stainless steel basket for equipment cleaning.
- ✓ A glass for cleaning of little instruments and its support.
- ✓ A lid that eliminates the noise while working and prevents dirt coming into the cleaner.
- ✓ Cleaner tweezers (Ref. 070520)

Height	220 mm
Width	425 mm
Depth	270 mm
Weight	4,5 kg
Tank capacity	4 l
Voltage	230 V - 50/60 Hz
Transducer power	100 W
Frecuency	40 khz

(L) 420 x (W) 360 x (H) 300 mm | 9 kg



Ref. 100152 - 5,8 L

Inner dimensions:
Height 155 mm
Width 298 mm
Depth 152 mm

Ref. 100151 - 3,4 L

Inner dimensions:
Height 100 mm
Width 298 mm
Depth 152 mm

Ref. 100150 - 2,8 L

Inner dimensions
Height 100 mm
Width 238 mm
Depth 140 mm

Ultrasonic bath: 5,8l - 3,4l - 2,8l

Made of stainless steel: high resistance to wear and tear and long service life. Easy-to-use digital time control. It has a drain on the side. Time from 1 to 99 minutes. Temperature from 1 to 80°C. Suitable for cleaning 3D printed models with resin.



Reference	R-100152	R-100151	R-100150
Height	290 mm	290 mm	245 mm
Width	365 mm	365 mm	300 mm
Depth	180 mm	165 mm	165 mm
Weight	5,400 kg	5,2 kg	3,980 kg
Tank capacity	5,8 L	3,4 l	2,8 L
Tension	230 V 50/60 Hz	220-240 V, 50/60 Hz	230 V 50/60 Hz
Power Ultr.	180 W	180 W	120 W
Heating power	200 W	200 W	100 W
Frecuency	40 Khz	40 kHz	40 Khz

Cleaning and disinfection

Extraoral suction system

Ref. 100347

Position the mouthpiece next to the patient's mouth and avoid the spread of all types of aerosols, particles in suspension, vapors and odors in the clinic. Greatly reduce the risk of cross contamination between professionals and patients and protect respiratory health.

- ✓ Two power levels.
- ✓ 2 filter cartridges composed of:
 - Coarse filter for small particles.
 - Active carbon filter.
 - HEPA H13 filter.
- ✓ UV light.
- ✓ Plasma generator.
- ✓ Ergonomic shape and easy to move to have a good working position.
- ✓ Semi-rigid suction arm, self-supporting and easy to clean.
- ✓ Easy cleaning housing.
- ✓ Hours of use counter. *

Voltage	220 V, 50 Hz
Power	1000 W.
2 power levels	(50-100 %)
Noise	60 dB.
Flow	4000 L/min.
Máx. vacuum	-25 kPa.
UV lamp	254 nm
Plasma generator	5000 V
HEPA grade	H13 (0.3 µm, > 99.97 %).
Arm length	1.3 m.



***Observation for maintenance:**
It is recommended to replace the filter element after one year of use.

Power supply	230 V / 50 Hz
Power	1.000 W
Tank	5 l.
Spray volume	470 ml/minutes
Droplet size	Adjustable: 100 - 150 µm
Effective range	Adjustable: 3-8 m.



Purifier / Nebulizer

Ref. 100520

High efficiency purifier-nebulizer device, capable of atomizing the disinfectant solution creating a very fine mist that penetrates every corner.

- ✓ The droplet size can be adjusted, as well as the jet launch distance.
- ✓ Using the appropriate products, it is possible to efficiently disinfect large surfaces in just a few minutes (up to 100 m² in just 10 minutes).

- ✓ Its 5 l tank. allows spraying many square meters on a single charge.
- ✓ The reduced weight of the nebulizer -about 3 kg-, together with the practical banner tape that is included, allows a comfortable, efficient and safe use.
- ✓ The device is essential to disinfect medium-sized surfaces.



Portable dental unit

Ref. 100650

Designed to bring the operational capacity of a dental office complete dental care to any location or complex situation.

Its compact and highly integrated design contains all the essential elements in dental practice.

Power supply	220V - 240V AC / 50 Hz
Nominal power	600 W
Airflow out	118 L/min a 0 bar
Handpiece pressure	0,2 - 0,25 MPa (2,0 - 2,5 bar)
Aspiration pressure	0,06 - 0,1 MPa (0,6 - 1,0 bar)
Noise level	40 - 45 dB

Cleaning and disinfection

Air purifiers

Family of air purifiers of three different sizes. They act on the stale air in any room, eliminating dust, dirt and germs. After the treatment, the air is much cleaner and more breathable.

The purification process begins with a battery of filters of three different levels: bacterial HEPA, activated carbon and photocatalytic. Subsequently, a highperformance UV lamp continues the cleaning work.

Reference	100510	100512	100514
Power	36 W	45 W	72 W
Noise	≤40 dB(A)	≤45 dB(A)	≤45 dB(A)
Purificated volume	150 m³/h	280 m³/h	560 m³/h
Recommended purification area	25 m²	30-35 m²	60-70 m²
3-level filtration (HEPA, activated carbon, catalytic)	Yes	Yes	Yes
UV cleaning + Ionization	Yes	Yes	Yes
UV lamp frequency	254 nm	254 nm	254 nm
Ozone generator	Yes	Yes	—
Remote control	Yes	Yes	Yes
Indicators of: temperature, humidity, air quality, filter change	No	Yes	Yes

Purifier 25 m²

Ref. 100510



303 mm
170 mm
329 mm
4 kg

Purifier 30-35 m²

Ref. 100512



536 mm
206 mm
360 mm
6,3 kg

Purifier 60-70 m²

Ref. 100514



570 mm
265 mm
370 mm
8,5 kg



Ref. 100500

Ideal for small objects such as mobile phones, masks, key rings, coins, jewelry, etc. Highly recommended for personal use. With drying function through heating resistance.

120 x 296 x 195 mm | 1 kg



Ref. 100502

Medium-sized. Its camera allows to process the same type of articles as the previous model, but in greater quantity; or if you prefer, larger objects. Ideal for clinics and laboratories.

300 x 290 x 200 mm | 4,5 kg



440 mm
365 mm
405 mm
6 kg

Hygienic Storage Cabinet **Ref. 100505**

Continuous protection for instruments and work materials in clinics and laboratories.

- ✓ Ultraviolet (UV) light technology.
- ✓ Dual technology (UV + HEAT).
- ✓ Automated cycle control.
- ✓ Integrated drying function.
- ✓ Professional construction.

Interior volume	38 L
Chamber material	Acero inoxidable
Electrical power	460 W
Power supply	230 V, 50 Hz

UV cabinets

Reference	100500	100502
Voltage	110~240 V, 50 Hz	220 V, 50 Hz
Power	150 W	4 W
Timer	10 min.	0-60 min/continuous



Prime 221 - intramatic connection

Ref. 100237

Dimension	90 x 191 x 130 mm	Power	120 W
Weight	2,8 kg	Voltage	0~36 V
Max. speed	50000 rpm	Weight handpiece	205 g
Torque	5.8 N-cm		




Sense A3 - intramatic connection

Ref. 100239

Dimension	88 x 147 x 120 mm	Power	70 W
Weight	1,15 kg	Voltage	0~33 V
Max. speed	35000 rpm	Weight handpiece	215 g
Torque	3.0 N-cm		

Sirocco S

Ref. 080534



Induction suction "S" + LED tube

Ref. 080534AT



Basic configuration of the Sirocco Model S table-top suction unit, supplied with a flexible and extendable hose up to 700 mm for directing the suction to the work area. It allows small and delicate polishing and reworking operations.

- ✓ Easy to orient the suction hood.
- ✓ Built-in LED lighting.

Tube 700 mm • Ø 45 mm • Weight 3.8 kg

Induction suction "S" + laboratory review box

Ref. 080534AL



Sirocco S suction unit and a dental laboratory finishing box that is fully compatible with the benchtop unit, forming a very compact and balanced unit.

- ✓ Spacious and robust design.
- ✓ LED lamp for perfect illumination.
- ✓ Comfortable and practical armrests.

Height 365 mm • Width 345 mm
Depth 400 mm • Weight 5.5 kg

Alginate mixer

Ref. 080515

Useful tool when making alginate paste in the dental clinic. Its operation is easy: just place powder and water in the bowl and select the rotation speed (fast/slow) by pressing the button. While rotating, the user can mix the paste with a spatula till it reaches the appropriate consistency and texture.

The user can control visually the process at every stage, and this guarantees that the paste made is in optimal conditions for application. In order to make cleaning easier the bowl is detachable. It is factory supplied with a bowl and a special plastic spatula.

Dimensions	
250 x 175 x 230 mm (with bowl)	
Weight	3 kg
Power	60 W

Ref. 080515-02	
Bowl	750 cc
Diameter	130 mm
Height	110 mm
Weight	0.2 kg



CE Medical Device MDR 2017/745

(L) 300 x (W) 310 x (H) 195 mm | 3.9 kg

Automatic mixing alginate

Ref. 100191

Just fill the bowl with alginate and water, put it into the machine, press the button, and in a few seconds you will achieve a completely homogeneous mixture, ready for use. Silent, small, comfortable and above all economical and reliable.

Factory supplied with a bowl set, a plastic spatula and a practical water dispenser. It can also be used as plaster mixer with astonishing results.

The machine is factory supplied with four bowls, a practical water dispenser, a mixing stick, a spatula and a holder.



MAXIMUN MIXTURE
80 grs. Alginate + water

Dimensions	265 x 205 x 305 mm
Weight	17 kg
Voltage	230 V 50/60 Hz
Power	350 W

(L) 300 x (W) 420 x (H) 430 mm | 18,5 kg



Tight tank for alginate

Ref. 100189



Dimensions	220 x 150 x 190 mm
Weight	0,467 kg
Capacity	1 kg

(L) 220 x (W) 155 x (H) 200 mm
0,62 kg

Silicone mixer

Ref. 100360

Impression material mixer simple, fast, reliable and safe. With small size to be located anywhere in the laboratory. The cartridge insertion/removal operation is easily performed within a few seconds.

All this means reduction of time, material and the improvement of the quality of the works.


Compatible with almost cartridges type 380 ml present in the market. Suitable for 99% of silicones.

The mixer has three adjustable pouring speeds to cover the needs of each job. Its amazing sturdiness together with the power of its motors generate a totally fluid mass, with optimum viscosity and totally free of bubbles without any effort.

Environment temperature	5-40°C
Supply	240 V 50 Hz
Power	150 W
Push	450 kg ±10



290 mm
170 mm
250+150 mm | 7,1 kg



Vacuum forming machine

Ref. 080542

It is cheap, easy to use and has an attractive design. It has a powerful aspiration motor with automatic connection when the sheet holder goes down. The carbon heater allows a progressive and uniform heating of the material that is going to be formed. It can be used with square sheets. Is supplied with a plane model holder and with a basket that can be filled with grit (a can is included with the machine) so forming depth is reduced.

(L) 197 x (W) 267 x (H) 337 mm | 5,7 kg



Pressure former PlastPress

Ref. 080544

The most compact positive pressure machine on the market. Greater accuracy in the forming process and lower noise compared to vacuum machines. Practical and versatile, it accepts round and square sheets up to 6 mm thickness. It must be connected directly to a compressor (between 3 and 8 bar).

(L) 310 x (W) 270 x (H) 310 mm | 8,5 kg



Pressure former

Ref. 080543

The strong pressurization allows precise molding, regardless of the cut of the model, as well as various thicknesses (0.5 to 4.0 mm). Square and round sheets (127, 125, 120 mm) can be used.

Dimensions	260 x 410 x 545 mm	Power	450 W
Weight	15 kg	Pressure	4-6 bar



Drying oven for thermoforming sheets

Ref. 080118

A fan circulates the hot air producing an uniform drying.

- ✓ Draws out moisture and keeps dry thermoforming sheets, to avoid the formation of pores.
- ✓ Drying of plaster models, instrument, etc.
- ✓ 6 stainless steel trays.
- ✓ Stainless steel housing.
- ✓ 24 hour timer.

Dimensions	450 x 340 x 310 mm	Power	550-650 W
Weight	8,3 kg	Temperature	35-75°C

Thermoforming sheets

MESTRA

DOUBLE BAG:

- Metallized outer bag, vacuum-sealed, with self-seal.
- Inner bag that contains the sheets.

MATERIAL	MM	QTY	REFERENCE
CLEAR Duras Hart Hard Rígidas Rígides Rigide	0,8	10	080801-10
		25	080801-25
	1,0	10	080802-10
		25	080802-25
	1,5	10	080803-10
		25	080803-25
	2,0	10	080804-10
		25	080804-25
	1,0	10	080812-10
	1,5	10	080813-10
MOUTHGUARD Blandas Weich Soft Flexíveis Souples Morbide	2,0	10	080814-10
	3,0	5	080815-05

FOR ALL TYPES OF USES:

- Orthodontic aligners
- Orthodontic appliances.
- Dental aligners.
- Space maintainers.
- Occlusal splints (bruxism).
- Registration plates for surgical guide.
- Prosthetic post-treatment.
- Registration plates.
- Deprogrammers.
- Base plates for complete dentures.
- Individual trays.
- Surgical guides.
- Apnoea devices.
- Trays for whitening.
- Fluoride trays.
- Sports mouthguards.

Thermoforming sheets



Standard square and round thermoforming sheets. Made of high quality plastic with CE marking certification as medical device. There is no soot, nor fumes during the heating. As they do not contain recycled materials, they do not tear while subjected to vacuum.

Ref.	Type	Hardness	Appearance	Size	Thickness	Qty.
080574	Crystal	Hard	Transparent	126x126	0.3 mm	10
080575	Crystal	Hard	Transparent	126x126	0.75 mm	5
080576	Crystal	Hard	Transparent	126x126	1.0 mm	5
080577	Crystal	Hard	Transparent	126x126	1.5 mm	5
080578	Crystal	Hard	Transparent	126x126	2.0 mm	5
080579	Milky	Hard	Translucent	126x126	0.6 mm	10
080580	Soft	Soft	Translucent	131x131	1.0 mm	10
080581	Soft	Soft	Translucent	131x131	1.5 mm	10
080582	Soft	Soft	Translucent	131x131	2.0 mm	10
080583	Soft	Soft	Translucent	131x131	3.0 mm	5
080584	Soft	Soft	Translucent	131x131	4.0 mm	5

Ref.	Type	Hardness	Appearance	Size	Thickness	Qty.
080574R	Crystal	Hard	Transparent	Ø125	0.3 mm	10
080575R	Crystal	Hard	Transparent	Ø125	0.75 mm	5
080576R	Crystal	Hard	Transparent	Ø125	1.0 mm	5
080577R	Crystal	Hard	Transparent	Ø125	1.5 mm	5
080578R	Crystal	Hard	Transparent	Ø125	2.0 mm	5
080579R	Milky	Hard	Translucent	Ø125	0.6 mm	10
080580R	Soft	Soft	Translucent	Ø125	1.0 mm	10
080581R	Soft	Soft	Translucent	Ø125	1.5 mm	10
080582R	Soft	Soft	Translucent	Ø125	2.0 mm	10
080583R	Soft	Soft	Translucent	Ø125	3.0 mm	5
080584R	Soft	Soft	Translucent	Ø125	4.0 mm	5




CE Medical médical

CLOSING METHODS:

In addition to using the pouch's zipper closure, it is advisable to double-fold the edge and secure it with clamps or clips to ensure a hermetic seal.




CE Medical médical



Endodontics motor

(wireless) **Ref. 100350**

- ✓ Optimized design for user interface.
- ✓ Ergonomic and compact body.
- ✓ 5 memory programs.
- ✓ From 100 to 600 r.p.m.
- ✓ Auto and Autostop inverse functions.



Maximum	100-600 rpm
Torque	1-4 N·cm
Voltage	230 W
Power supply	Pilas

(L) 240 x (W) 170 x (H) 80 mm
0,7 kg



Apex locator


Ref. 100330

- ✓ Automatic calibration, LCD display clear and accurate image.
- ✓ Different colours.
- ✓ Sterilizable accessories.
- ✓ Rechargeable battery.

Dimensions 80 x 138 x 110 mm

Weight 0,4 kg

(L) 230 x (W) 155 x (H) 112 mm | 0,8 kg



Gutta-Percha Cutter

Ref. 100349

- ✓ The tip is flexible design to let you cut smoothly.
- ✓ Press the power button and the tip can heat up to 500 °C immediately.
- ✓ Release the power button and the tip can cool down fast.
- ✓ Put tip cap back after using.
- ✓ Put GP Cut back on stand base.




Voltage	2 AA batteries (DC 3V)
Temp. Reach	>500 °C
Cool down period	Less than 1 second
Stand base size	3 cm x 2,5 cm
3 tips	(L: 2,0 mm - M: 1,2 mm - S: 0,5 mm)

(L) 180 x (W) 70 x (H) 35 mm
0,12 kg



2 batteries
AA (3V cc)

- ✓ Wide lighting range (from 420 to 472 nm).
- ✓ High power (up to 2000 mW/cm²).
- ✓ Easily adjustable in height and orientation.
- ✓ Digital display.
- ✓ Remote control.
- ✓ It is factory supplied complete, with 20 disposable mouth splints, 2 protection goggles, and a lip retractor.



Whitening lamp

Ref. 100340

- ✓ Wide spectrum of illumination.
- ✓ 5 Blue lights, 2 red lights.
- ✓ 3 different types of LED light sources: blue, red, (red + blue).
- ✓ 8 purple lights with sterilization effect.
- ✓ Made of metallic material, easy to assemble.
- ✓ Adjustable bleaching time.




Whitening lamp

Ref. 100341



Mixing frequency	5000 rpm
Timer range	4 ~ 60 s
Motor power	40 W
Power supply	230 V, 50 Hz



Capsule mixer

Ref. 100625

- ✓ High speed.
- ✓ Universal compatibility.
- ✓ Precision digital control.

- ✓ Ergonomic and safe design.
- ✓ Stability and quiet operation.
- ✓ Durability.



Intraoral camera

(with autofocus system)

Ref. 100336

- ✓ 21.5-inch LCD touch screen.
- ✓ High definition camera (4K).
- ✓ Autofocus.
- ✓ Image and video recording on USB memory.
- ✓ QR code for downloading images via WIFI.
- ✓ Support bar.
- ✓ Android system.
- ✓ Display of 1 or 4 images on screen.




Composite heater

Ref. 080077

Device to keep composite syringes warm.

Temperature	40-50 y 70 °C
Ø of holes	16 mm
Power	100 a 240 V



Screen: 2,52 kg
Camera: 0,1 kg



Made in Spain

Chair for workstation

Wheels with safety brake (EN 1335-2:2009)

Ref. 100300



(Kit) Finishing table


Ref. 100301 K2

- ✓ Metal table, wood worktop.
- ✓ Chair for workstation.
- ✓ Polishing lathe + 2 trays.
- ✓ Micromotor.
- ✓ LED lamp.

GIFT:
→ Magnifying glass.
→ Magnet bur holder.

(L) 950 x (W) 610 x (H) 100 mm | 63 kg

Made in Spain



- 1 100301-04 Aspiration peg.
- 2 100301-05 Peg support.
- 3 100301-02 Metal tube.
- 4 080530-04 Flexible tube (1 m.)




Finishing table

Ref. 100301

Metal work table with wooden worktop. Assembly is simple.

They have all the elements and accessories for work in dental environments: 2 sockets, shelf, drawer, etc.

(L) 900 x (W) 600 x (H) 660 mm | 46 kg



Work chairs

Work chairs specially designed for the dental laboratory.

Comfortable, ergonomic and easy to clean. Allow 1 seat height adjustment. The material used is spongy polyurethane in black.



	Ref. 100297	Ref. 100298
Max. seat height	660 mm	555 mm
Min. seat height	460 mm	420 mm
Seat depth	450 mm	450 mm
Seat width	460 mm	460 mm
Weight	9 kg	8 kg

Wheels with safety brake (EN 1335)



Work chair

Ref. 100299

4 different regulations allowed: seat height and inclination and backrest height and angle of inclination.

Max. height	555 mm
Min. height	420 mm
Depth	450 mm
Width	460 mm
Weight	10.3 kg

Wheels with safety brake (EN 1335)



Work chair

Ref. 100300

4 different regulations allowed: seat height and inclination and backrest height and angle of inclination.

Wheels with safety brake (EN 1335-2)

Max. seat height	515 mm
Min. seat height	390 mm
Seat depth	420 mm
Seat width	420 mm
Weight	9 kg

(L) 610 x (W) 610 x (H) 325 mm | 10,5 kg



Poni chair without backrest

Ref. 100600



Wheels with safety brake (EN 1335)

Max. seat height	730 mm
Min. seat height	530 mm
Seat depth	320 mm
Seat width	360 mm
Weight	7 kg




Poni chair with backrest

Ref. 100602

- ✓ Height Adjustable.
- ✓ High quality PU leather.
- ✓ High resilience foam.
- ✓ Four-stage inlet pneumatic bar. Safe + Smooth. TUV and SGS Authentication.
- ✓ Hydraulic oil seal. Multiple protection, no oil leakage, no air leakage. Long service life.
- ✓ Oil filter cleanliness. Built-in oil scraping device.



Wheels with safety brake (EN 1335)



2 wheels trolley

Ref. 100451


- ✓ Made of stainless steel.
- ✓ Removable tray (495 x 345 mm)



5 wheels trolley

Ref. 100452

- ✓ Made of stainless steel.
- ✓ Removable tray (495 x 345 mm)



Scanner trolley

Ref. 100450

- ✓ Material: antibacterial plastic, aluminum, and steel.



Intraoral scanner trolley with CPU support

Ref. 100449

- ✓ Material: antibacterial plastic, aluminum, and steel.




Electric mobile trolley

Ref. 100305

Made of steel sheet and enameled in white.

- ✓ 3 trays of 40 x 30 cm.
- ✓ 2 sockets to power the equipment.
- ✓ 4 wheels with smooth and balanced movement.



(L) 510 x (W) 420 x (H) 840 mm | 15 kg




Electric mobile trolley

Ref. 100306

Made of steel sheet and enameled in white.

- ✓ 3 glass trays.
- ✓ 2 sockets to power the equipment.
- ✓ 4 wheels with smooth and balanced movement.



(L) 575 x (W) 560 x (H) 845 mm | 25,5 kg




Electric mobile trolley

Ref. 100307

Made of steel sheet and enameled in white.

- ✓ 3 trays of 290 x 396 mm
- ✓ 3 sockets to power the equipment.
- ✓ 4 wheels with smooth and balanced movement.
- ✓ Weight supported per tray: 5 kg



(L) 430 x (W) 520 x (H) 860 mm | 12 kg




Electric mobile trolley

Ref. 100308

Made of steel sheet and enameled in white.

- ✓ 3 trays of 290 x 396 mm.
- ✓ 3 sockets.
- ✓ 1 drawer.
- ✓ 4 wheels with smooth and balanced movement.
- ✓ Weight supported per tray: 5 kg



(L) 430 x (W) 520 x (H) 860 mm | 16 kg




Electric mobile trolley

Ref. 100309

NEW

Electric mobile trolley with wheels for the dental clinic. Made of resin and steel sheet enameled in white.

- ✓ 3 sockets to power the equipment.
- ✓ 3 drawers.
- ✓ 1 power sockets.
- ✓ White colour.
- ✓ Smooth glide.



(L) 480 x (W) 540 x (H) 840 mm | 21 kg




Assemblable mobile trolley

Ref. 100307M

NEW

Made of steel sheet and enameled in white.

- ✓ 3 trays of 290 x 396 mm
- ✓ 3 sockets to power the equipment.
- ✓ 4 wheels with smooth and balanced movement.
- ✓ Weight supported per tray: 5 kg



REDUCED PACKAGING:
(L) 460 x (W) 780 x (H) 120 mm | 11,5 kg



Trolley support

Ref. 100450-01

NEW

Waste trolley

Ref. 100440

Ref. 100442

Two mobile containers designed for the safe transport of infectious medical waste. Thanks to their simple and functional design, they fit into any type of environment. Waste trolley

Available for 40 L or 60 L.

Ref. 100440 - 40 L	
Height 568 mm	Width 296 mm
Depth 404 mm	Weight 11 kg

Ref. 100442 - 60 L	
Height 618 mm	Width 346 mm
Depth 454 mm	Weight 14 kg



Water bath

Ref. 080078

This machine is extremely useful when working with applications that need a water bath that keeps the water temperature accurately (polymerisation of small elements, wax sheets for biting registers, etc). It has a 1 litre stainless steel tank and it is easy to clean. The body is made in high density polyurethane, resistant to chemical agents.

Bath temperature control is made by an electronic thermostat easy to use and accurate when keeping the programmed temperature.

Temperature range is from 30°C to 75°C. A powerful heater allows reaching the set temperature in a few minutes.

Height	146 mm	Voltage	230 V
Width	230 mm	Temperature range	30-75 °C
Depth	195 mm	Power	310 W
Weight	1.7 kg	Capacity	1 l.

(L) 300 x (W) 300 x (H) 200 mm | 2 kg



Dry piston compressors

NOT SUITABLE FOR MILLING MACHINES.

Family of four compressors consisting of two sizes of 80 and 160 litres/minute to which a dryer with automatic drain can be coupled.

They are supplied with two pressure gauges (boiler and service pressure), outlet pressure regulator and complete set of valves for draining. The service pressure can reach 8 bars.

Extremely quiet operation (60 dB), easy maintenance (dry piston) and excellent reliability.

Ref. 110300 - Compressor: 80 l/min
Ref. 110310 - Compressor: 80 l/min (with dryer)
Ref. 110330 - Compressor: 160 l/min
Ref. 110340 - Compressor: 160 l/min (with dryer)




How do I know if my compressor complies with the European safety regulations?

It is easy to find dental compressors made in the East that do not comply with current European safety regulations.

For this reason, it is advisable to check its features and security stamp: it is a plate welded to the tank similar to the attached picture. Next to the CE mark, you can see the number of the Notified Body that has validated its adaptation to European standards.

1. Ref. 110340 is also available with HEPA filter (H14) Ref. 110340H
2. The activated carbon filter Ref. 110300-53 is available as an accessory for all compressors.



Reference	110300	110310	110330	110340
Power	750 W	750 W	1500 W	1500 W
Effective airflow (4 bars)	80 l/min	80 l/min	160 l/min	160 l/min
Noise	53 dB	53 dB	54 dB	54 dB
Pressure	5-8 bar	5-8 bar	5-8 bar	5-8 bar
Tank capacity	30 l	30 l	50 l	50 l
Heads	1	1	2	2
Depth	390 mm	540 mm	370 mm	480 mm
Width	390 mm	390 mm	680 mm	680 mm
Height	650 mm	650 mm	750 mm	750 mm
Weight	30 kg	37 kg	58 kg	66 kg

CERTIFICATE by TÜV SÜD Type certificate N° 12202200491258001
Directiva 2014/29/UE (pressure vessels)

Wet ring suction unit

Ref. 100345

Designed for a simple and reliable operation, the unit gives a quick response: it detects the start-up signal and reaches the optimum vacuum level in a record time. Its operation is extremely silent and it does not need any maintenance, since this technology automatically separates the air and the water with debris.


The suction power is always kept stable, with no noticeable ups and downs or failures. The unit is useful for intensive use workplaces, or for two workplaces in case of moderated use.

(L) 470 x (W) 430 x (H) 650 mm | 25 kg

Dimensions	360 x 405 x 585 mm
Weight	27 kg
RPM	0-4000
Air flow	300 l/min
Power	750 W
Noise	58 dB




Lightning



LED lamp (articulated arm)

Ref. 100288


- ✓ 84 large LEDs.
- ✓ Color temperature 6000 - 7000 K
- ✓ Power consumption 2 to 21W
- ✓ Touch-adjustable light intensity.
- ✓ Free of Flicker effect.
- ✓ With clamp for fastening to the table.
- ✓ Highly resistant articulated arm.



LED lamp (articulated arm)

Ref. 100290

- ✓ 117 small LEDs.
- ✓ Color temperature 600 - 6000 K
- ✓ Power consumption 6 to 24 W
- ✓ Free of Flicker effect.
- ✓ With clamp for fastening to the table.



LED lamp (articulated arm)

Ref. 100289

- ✓ 90 large LEDs.
- ✓ Color temperature 6000 - 7000 K
- ✓ Touch-adjustable light intensity.
- ✓ With clamp for fastening to the table.




Magnifying glass with LED lamp

Ref. 100287 (3 diopter lens)

Ref. 100287-5x (5 diopter lens)


- ✓ 90 white LEDs.
- ✓ 3 - 5 diopter lens.
- ✓ Ø 120 mm lens.
- ✓ Clamp to fix it to the table.



LED lamp (flexible arm)

Ref. 100289V

- ✓ 180 LEDs
- ✓ Temperature color 6500 - 200 K
- ✓ Energy consumed 30 W
- ✓ 3 types of light: cold, neutral and warm.
- ✓ Tactile adjustable light intensity.
- ✓ Flicker-free.
- ✓ With table clamp.



Accessories for lamps and magnifying glass

Ref. 100287 - 100287-5X - 100288 - 100290

Alphabetical index

Description	Page
ADERER PLIERS.....	4
AIR PURIFIER	27
ALCOHOL LAMP.....	11
ALGINATE BOWL.....	3
ALGINATE CONTAINER.....	29
ALGINATE MIXER.....	29
ALGINATE SPATULA STAINLESS STEEL.....	3
AMALGAM APPLIER.....	5
ANGLE PLIERS.....	4
APEX LOCATOR DPEX III.....	31
ASPIRATION, SIROCCO.....	28
ARC FOR COFFERDAM HOLDING.....	9
AUTOCLAVE.....	17
AUTOSEALING BAGS FOR SEALING MACHINE.....	24
BAGS FOR SEALING MACHINE, AUTOSEALING.....	24
BARRIER FILM.....	24
BASE MODELFORMER.....	3
BIOART MAGNIFYING GLASSES.....	7
BIOLOGICAL INDICATOR INCUBATOR.....	19
BOBBIN FOR SEALING MACHINE.....	24
BOX FOR DENTURES.....	13
BOX FOR ENDODONTIA.....	13
BOX FOR ORTHODONTIC MODELS.....	14
BOX FOR ORTHODONTIC PLATES.....	13
BOX FOR PROTHESIS.....	14
BOXES.....	13
BOXES FOR MODELS.....	14
BOXES FOR TRANSPORT.....	14
BRUSHES FOR PROPHYLAXIS.....	10
BUR CLEANER TWEEZERS.....	9
BUR STAND.....	10
CABINETS, UV.....	27
CALIPER.....	6
CALYPSO.....	22
CAMERA, INTRAORAL.....	31
CARVERS PK THOMAS.....	4
CAPSULE MIXER.....	32
CEMENT SPATULA.....	4
CHAIR FOR WORK.....	33
CHEEK RETRACTOR.....	9
CLAMPS HOLDER.....	6
CLIPS FOR DENTAL FILM.....	8
COFFERDAM HOLDING ARC.....	9
COMPOSITE HEATER.....	31
COMPRESSORS.....	36
CONTAINER FOR ALGINATE.....	30
CONTAINER FOR MODELS.....	13
CORD.....	13
COTTON CONTAINER.....	11
COTTON ROLL DISPENSER.....	11
CROWN REMOVING HAMMER.....	5
CUPS FOR PROPHYLAXIS.....	9
CUTTER FOR GUTTA PERCHA.....	31
CURING LAMP.....	21
DAM PERFORATOR.....	4
DAPPEN GLASSES.....	3
DECALCIFIER.....	19
DENTAL LOUPES.....	7
DENTAL FILM CLIPS.....	8
DEVELOPER, MANUAL.....	6
DISINFECTANT, NEBULIZER.....	27
DISPENSER FOR COTTON ROLLS.....	11
DISPENSER OF GLASSES.....	11
DISTILLED WATER, TANK.....	19
DISTILLER - DECALCIFIER.....	19
DOSE MEASURING SPOON.....	8
DOUBLE GUIDE OCCLUSOR.....	3
DRYING OVEN.....	30

Description	Page
ENDODONTIC MOTOR.....	31
ENDODONTIA BOX.....	13
ENDODONTIA, RULE FOR MEASURING FILES.....	13
EXTRAORAL SUCTION UNIT.....	25
FACE BOW.....	3
FACE PROTECTION.....	7
FILM CLIPS.....	8
FINISHING TABLE.....	33
FLUORIDE FOAM.....	13
FORMING SHEETS.....	29
FOX PLANE.....	9
GERMICIDE BOX (FOR BURS).....	10
GLASSES, DAPPEN.....	3
GLASSES, DISPENSER.....	11
GUN FOR SANDBLASTING.....	19
GUTTA-PERCHA CUTTER.....	31
HAMMER FOR CROWN REMOVING.....	5
HANDLE FOR MIRRORS.....	5
HANDLE FOR MORSE POINT.....	5
HEAD MAGNIFIER.....	7
HEATER FOR COMPOSITES.....	31
HOLDER FOR INSTRUMENTS.....	3
HOLDING CLAMPS FOR WELDING.....	6
HYGIENIC STORAGE CABINET.....	27
IMPRESSION TRAYS.....	13
INCUBATOR FOR BIOLOGICAL INDICATOR.....	19
INSTRUMENT TRAY.....	16
INTRAORAL CAMERA.....	31
IVORY TYPE DAM PERFORATOR.....	4
LAMP FOR ALCOHOL.....	11
LAMP FOR TEETH WHITENING.....	31
LED LAMP.....	36
LIQUID DOSIFIER.....	11
LOCATOR, APEX.....	31
MAGNETIC INSTRUMENTS RACK SUPPORT.....	3
MAGNETIC STAND.....	10
MAGNIFIER FOR HEAD.....	7
MAGNIFIER GLASSES.....	7
MAGNIFYING GLASS.....	36
MAGNIFYING GLASSES, BIOART.....	7
MATHIEU PLIERS.....	4
MATRIX RETAINERS.....	5
MICROAPPLICATORS.....	9
MICROMOTORS INTRAMATIC.....	28
MINI SANDBLASTING.....	19
MIRROR HANDLE.....	5
MIRRORS.....	5
MIXER FOR SILICONE.....	30
MIXING MORTAR.....	5
MIXING TIPS.....	12
MODEL CONTAINER.....	13
MODEL FORMER.....	3
MODEL TRAY.....	15
MORSE POINT HANDLE.....	5
MORTAR FOR MIXING.....	5
MOTOR FOR ENDODONTIC.....	31
NAPKIN HOLDER.....	9
NEBULIZER FOR DISINFECTANT.....	27
ORTHODONTIC PLATES BOX.....	13
OVEN FOR DRYING.....	28
PK THOMAS CARVERS.....	4
PLASTER SPATULA.....	4
PLASTPRESS.....	30
PLATES FOR FORMING.....	3
PLIERS RACK SUPPORT.....	4
POLIMERIZATION LAMP LED.C.....	21
PRECISION SCALE.....	8
PRESSURE FORMER.....	28

Alphabetical index

<u>Description</u>	<u>Page</u>
PRESSURE FORMER - PLASTPRESS	28
PROTHESIS BOX	13
PROPHYLAXIS, BRUSHES	10
PROPHYLAXIS, CUPS	10
PURIFIER FOR AIR	26
RAIL FOR MODEL TRAY	15
REEL FOR SEALING MACHINE	24
RETRACTOR FOR CHEEKS	9
RULE FOR MEASURING FILES	16
SAFETY GOGGLES	7
SANITIZER, UV	27
SCHWARZ PLIERS	4
SEALING MACHINE	24
SELF-SEALING STERILISATION POUCH	24
SILICONE MIXER	29
SILICONE PISTOL	12
SIROCCO, ASPIRATION	28
SPATULA FOR ALGINATE	3
SPATULA FOR CEMENT	4
SPATULA FOR PLASTER	3
SPOON FOR DOSE MEASURING	8
STAND FOR BURS	10
STERILIZATION, TONGS	5
STERILIZATION, BOX FOR BURS	10
SUCTION UNIT, EXTRAORAL	25
SUCTION UNIT, WET RING	32
SUPPORT FOR PLIERS	4
SUPPORT FOR INSTRUMENTS, MAGNETIC RACK	3
SUPPORT RAIL FOR MODEL TRAY	15
TABLE FOR FINISHING	33
TANK FOR DISTILLED WATER	19
TEETH WHITENING LAMP	31
TERMOPLAST	30
THERMOFORMING SHEETS	31
TIPS FOR MIXING	12
TONGS FOR STERILIZATION	5
TOYS	10
TRANSPORT BOXES	14
TRAYS FOR IMPRESSION	12
TROLLEY	35
TROLLEY SUPPORT	35
TWEEZERS BUR-CLEANER	5
ULTRASONIC CLEANERS	22-23
UV RAY SAFETY SCREEN	7
UV SANITIZER CABINET	27
UV, PROTECTION GOGGLES	7
VACUUM FORMER	30
WALDSACHS PLIERS	4
WASTE TROLLEY	35
WATER BATH	36
WAX MODELING CARVERS PK THOMAS	4
WEINGART PLIERS	4
WELDING, HOLDING CLAMPS	6
WET RING SUCTION UNIT	32
WHITENING, TEETH LAMP	31
WORK CHAIR	33
X-RAY VIEWER	8



MESTRA

Talleres Mestraitua SL
(+34) 944 530 388
info@mestra.es - www.mestra.es
48150 SONDIKA (Spain)



Management
System
ISO 9001:2015

www.tuv.com
ID 1100018003